



PRECISIONARY

PRECISIONARY INSTRUMENTS:

ROTARY MICROTOME RF-800 For Histology & Pathology Solutions

Semy-Automated Rotary Microtome



Key Features

- 1** LCD touchscreen panel shows all cutting metrics.
- 2** Touchscreen control panel provides easy and friendly controls.
- 3** Specimen clamp is compatible with multiple cassette sizes.
- 4** Cutting stage has blade ejector and blade guard for increased safety. Can be used with high- or low-profile cutting blades.
- 5** Large volume waste tray is removable for easy cleaning.
- 6** Top of rotary microtome can be used as storage for consumables during cutting.
- 7** Smooth turning hand wheel can be locked anywhere in the 360° turn, allowing for improved ergonomics.

